

# Signia RIC product style and color guide

## Government Services

November 2020

**signia**  
Life sounds brilliant.

### Stylecto™ X

NEW



Made for  
iPhone | iPad | iPod



Rechargeable



Black/Silver   Cosmic Blue/Rose Gold   Snow White/Rose Gold   Black   White

### Pure® Charge&Go X



Made for  
iPhone | iPad | iPod



Rechargeable


### Pure 312T X



Made for  
iPhone | iPad | iPod



312 Battery Type



Beige   Granite   Grey   Black   Silver   Pearl White   Sandy Brown   Dark Champagne   Deep Brown   Rose Gold


Actual colors may vary.

### Click Sleeve 2.0


Products	Content	Item No.	
Type vented	Size		
	xS	10951152	
	S	10951146	
	M	10951148	
L		10951150	
	closed	xS	10951153
		S	10951147
		M	10951149
L		10951151	

### Fitting Ranges

S-Receiver Fitting Range  
45/108 - 45 dB Receiver




M-Receiver Fitting Range  
60/119 - 60 dB Receiver




### Click Domes

Products	Content	Item No.	
Type Open	Size		
	4 mm	10426016	
	8 mm	10426014	
10 mm		10426015	
	Semi-Open	8-12 mm	10670790
Closed		6 mm	10426017
		8 mm	10426012
		10 mm	10426013
Double Closed	8/10 mm	10426018	
	10/12 mm	10426019	


P-Receiver Fitting Range  
70/124 - 70 dB Receiver



Integrated HP-Receiver Fitting Range  
75/130 - 75 dB Receiver  
(Not available for Stylecto X)



## Common RIC Earmolds



		RIC	Embedded RIC <sup>1</sup>	RIC Lock	Embedded RIC Lock <sup>1</sup>	RIC Foil <sup>2</sup>	Foil RIC with Lock <sup>2</sup>	Canal RIC	Canal Lock RIC
	ROES Codes	EM11	EM28	EM12	EM29	EM15	EM16	EM17	EM18
Material	Acrylic (Hard)	•	•	•	•	•	•	•	•
	Silicone (Soft)	•	•	•	•	-	-	•	•
Vent Type	No Vent	•	•	•	•	•	•	•	•
	Semi-IROS	•	•	•	•	•	•	•	•
	Standard	•	•	•	•	•	•	•	•
Vent Size	1.0mm	•	•	•	•	•	•	•	•
	1.2mm	•	•	•	•	•	•	•	•
	1.6mm	•	•	•	•	•	•	•	•
	2.0mm	•	•	•	•	•	•	•	•
	3.0mm	•	•	•	•	•	•	•	•
	Largest Vent Possible	•	•	•	•	•	•	•	•

1) Embedded fixed receiver 2) Acrylic material only, applies to the Gloss

The **Bluetooth**® word mark and logos are owned by the **Bluetooth** SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners. Android is a trademark of Google Inc. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license. "Made for iPhone" means that an electronic accessory has been designed to connect specifically to iPhone and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPhone may affect wireless performance.

## Common Earmold Color Options

### Gloss\*



### Matte



### Opaque\*



\*Acrylic and Silicone for Black and White, Silicone only for other Opaque colors